

Product Datasheet

Rat Matrix Gla Protein (MGP) ELISA Kit (orb1203774)

Description	Rat Matrix Gla Protein (MGP) ELISA Kit
Reactivity	Rat
Range	93.7-6000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08494
Sensitivity	39.65pg/mL
Entrez	25333
Expiration Date	Please enquire.