

Product Datasheet

Rat Low Density Lipoprotein Receptor (LDLR) ELISA Kit (orb1203785)

Description	Rat Low Density Lipoprotein Receptor (LDLR) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P35952
Sensitivity	0.65ng/mL
Entrez	300438
Expiration Date	Please enquire.