

Product Datasheet

Rat Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit (orb1203797)

Description	Rat Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P0DMW1
Sensitivity	0.54ng/mL
Entrez	24472
Expiration Date	Please enquire.