

Product Datasheet

Rat Endothelin Receptor B (ETRB) ELISA Kit (orb1203831)

Description	Rat Endothelin Receptor B (ETRB) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P21451
Sensitivity	0.06ng/mL
Entrez	50672
Expiration Date	Please enquire.