

Product Datasheet

Rat Cysteinyl Leukotriene Receptor 2 (CYSLTR2) ELISA Kit (orb1203842)

Description	Rat Cysteinyl Leukotriene Receptor 2 (CYSLTR2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q924T9
Sensitivity	0.108ng/mL
Entrez	170926
Expiration Date	Please enquire.