

Product Datasheet

Rat Breast And Kidney Expressed Chemokine (BRAK) ELISA Kit (orb1203867)

Description	Rat Breast And Kidney Expressed Chemokine (BRAK) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8K453
Sensitivity	0.058ng/mL
Entrez	306748
Expiration Date	Please enquire.