

Product Datasheet

Rat A Disintegrin And Metalloproteinase With Thrombospondin 1 (ADAMTS1) ELISA Kit (orb1203900)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 1 (ADAMTS1) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q68EJ2
Sensitivity	30.5pg/mL
Entrez	79252
Expiration Date	Please enquire.