

Product Datasheet

Rat A Disintegrin And Metalloprotease 17 (ADAM17) ELISA Kit (orb1203902)

Description	Rat A Disintegrin And Metalloprotease 17 (ADAM17) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Z1K9
Sensitivity	0.112ng/mL
Entrez	57027
Expiration Date	Please enquire.