

Product Datasheet

Mouse VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit (orb1203928)

Description	Mouse VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q5UW37
Sensitivity	5.8pg/mL
Entrez	232983
Expiration Date	Please enquire.