

Product Datasheet

Mouse Uroporphyrinogen Decarboxylase (UROD) ELISA Kit (orb1203931)

Description	Mouse Uroporphyrinogen Decarboxylase (UROD) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P70697
Sensitivity	0.056ng/mL
Entrez	22275
Expiration Date	Please enquire.