

Product Datasheet

Mouse Uncoupling Protein 1, Mitochondrial (UCP1) ELISA Kit (orb1203934)

Description	Mouse Uncoupling Protein 1, Mitochondrial (UCP1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P12242
Sensitivity	0.061ng/mL
Entrez	22227
Expiration Date	Please enquire.