

Product Datasheet

Mouse Urocortin 2 (UCN2) ELISA Kit (orb1203936)

Description	Mouse Urocortin 2 (UCN2) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99ML8
Sensitivity	0.52ng/mL
Entrez	171530
Expiration Date	Please enquire.