


Product Datasheet

pCardiac Troponin I antibody (orb1204)

Description	pCardiac Troponin I antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human cardiac Troponin I around the phosphorylation site of Thr143 (RP(p-T)LR)
Target	phospho-cardiac Troponin I (Thr143)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P19429
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, IF=1:50-200, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow
cytometry of
K562.Primary
Antibody ...