

Product Datasheet

Mouse Sclerostin Domain Containing Protein 1 (SOSTDC1) ELISA Kit (orb1204001)

Description	Mouse Sclerostin Domain Containing Protein 1 (SOSTDC1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9CQN4
Sensitivity	0.068ng/mL
Entrez	66042
Expiration Date	Please enquire.