

Product Datasheet

Mouse Suppressors Of Cytokine Signaling 2 (SOCS2) ELISA Kit (orb1204003)

Description	Mouse Suppressors Of Cytokine Signaling 2 (SOCS2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	035717
Sensitivity	0.066ng/mL
Entrez	216233
Expiration Date	Please enquire.