

Product Datasheet

Mouse Suppressors Of Cytokine Signaling 1 (SOCS1) ELISA Kit (orb1204004)

Description	Mouse Suppressors Of Cytokine Signaling 1 (SOCS1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O35716
Sensitivity	0.14ng/mL
Entrez	12703
Expiration Date	Please enquire.