

Product Datasheet

Mouse Synaptosomal Associated Protein 25kDa (SNAP25) ELISA Kit (orb1204005)

Description	Mouse Synaptosomal Associated Protein 25kDa (SNAP25) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P60879
Sensitivity	0.14ng/mL
Entrez	20614
Expiration Date	Please enquire.