

Product Datasheet

Mouse Mothers Against Decapentaplegic Homolog 3 (Smad3) ELISA Kit (orb1204007)

Description	Mouse Mothers Against Decapentaplegic Homolog 3 (Smad3) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8BUN5
Sensitivity	0.136ng/mL
Entrez	17127
Expiration Date	Please enquire.