

Product Datasheet

Mouse Mothers Against Decapentaplegic Homolog 2 (Smad2) ELISA Kit (orb1204008)

Description	Mouse Mothers Against Decapentaplegic Homolog 2 (Smad2) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q62432
Sensitivity	24.6pg/mL
Entrez	17126
Expiration Date	Please enquire.