

Product Datasheet

Mouse Slit Homolog 3 (Slit3) ELISA Kit (orb1204009)

Description	Mouse Slit Homolog 3 (Slit3) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9WVB4
Sensitivity	14.5pg/mL
Entrez	20564
Expiration Date	Please enquire.