

Product Datasheet

Mouse Succinate Dehydrogenase Complex Subunit C (SDHC) ELISA Kit (orb1204035)

Description	Mouse Succinate Dehydrogenase Complex Subunit C (SDHC) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9CZB0
Sensitivity	0.071ng/mL
Entrez	66052
Expiration Date	Please enquire.