

Product Datasheet

Mouse Sphingosine 1 Phosphate Receptor 1 (S1PR1) ELISA Kit (orb1204043)

Description	Mouse Sphingosine 1 Phosphate Receptor 1 (S1PR1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O08530
Sensitivity	0.115ng/mL
Entrez	13609
Expiration Date	Please enquire.