

Product Datasheet

Mouse S100 Calcium Binding Protein A4 (S100A4) ELISA Kit (orb1204044)

Description	Mouse S100 Calcium Binding Protein A4 (S100A4) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07091
Sensitivity	14.2pg/mL
Entrez	20198
Expiration Date	Please enquire.