

## Product Datasheet

### Mouse C-Ros Oncogene 1, Receptor Tyrosine Kinase (ROS1) ELISA Kit (orb1204048)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| <b>Description</b>       | Mouse C-Ros Oncogene 1, Receptor Tyrosine Kinase (ROS1) ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                             |
| <b>Range</b>             | 78.1-5000pg/mL                                                    |
| <b>Note</b>              | For research use only                                             |
| <b>Application notes</b> | sandwich                                                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids.     |
| <b>Uniprot ID</b>        | <b>Q78DX7</b>                                                     |
| <b>Sensitivity</b>       | 33pg/mL                                                           |
| <b>Entrez</b>            | <b>19886</b>                                                      |
| <b>Expiration Date</b>   | Please enquire.                                                   |