

## Product Datasheet

### Mouse Replication Initiator 1 (REPIN1) ELISA Kit (orb1204054)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| <b>Description</b>       | Mouse Replication Initiator<br>1 (REPIN1) ELISA Kit |
| <b>Reactivity</b>        | Mouse                                               |
| <b>Range</b>             | 0.312-20ng/mL                                       |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | sandwich                                            |
| <b>Sample Types</b>      | tissue homogenates or<br>other biological fluids.   |
| <b>Uniprot ID</b>        | <b>Q5U4E2</b>                                       |
| <b>Sensitivity</b>       | 0.135ng/mL                                          |
| <b>Entrez</b>            | <b>58887</b>                                        |
| <b>Expiration Date</b>   | Please enquire.                                     |