

## Product Datasheet

### Mouse Regenerating Islet Derived Protein 3 Beta (REG3b) ELISA Kit (orb1204057)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Regenerating Islet<br>Derived Protein 3 Beta<br>(REG3b) ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                                   |
| <b>Range</b>             | 9.37-600pg/mL                                                           |
| <b>Note</b>              | For research use only                                                   |
| <b>Application notes</b> | sandwich                                                                |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids.                            |
| <b>Uniprot ID</b>        | <b>P35230</b>                                                           |
| <b>Sensitivity</b>       | 4.1pg/mL                                                                |
| <b>Entrez</b>            | <b>18489</b>                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                         |