

Product Datasheet

Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit (orb1204060)

Description	Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Reactivity	Guinea pig
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Uniprot ID	P97272
Sensitivity	0.059ng/mL
Entrez	100135495
Expiration Date	Please enquire.