

Product Datasheet

Mouse Peroxiredoxin 1 (PRDX1) ELISA Kit (orb1204074)

Description	Mouse Peroxiredoxin 1 (PRDX1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P35700
Sensitivity	0.103ng/mL
Entrez	18477
Expiration Date	Please enquire.