

Product Datasheet

Mouse Phosphatidylethanolamine Binding Protein 1 (PEBP1) ELISA Kit (orb1204093)

Description	Mouse Phosphatidylethanolamine Binding Protein 1 (PEBP1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P70296
Sensitivity	0.055ng/mL
Entrez	23980
Expiration Date	Please enquire.