

Product Datasheet

Mouse Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit (orb1204096)

Description	Mouse Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P08003
Sensitivity	0.112ng/mL
Entrez	12304
Expiration Date	Please enquire.