

Product Datasheet

Mouse Phosphodiesterase 1C, Calmodulin Dependent (PDE1C) ELISA Kit (orb1204098)

Description	Mouse Phosphodiesterase 1C, Calmodulin Dependent (PDE1C) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q64338
Sensitivity	0.053ng/mL
Entrez	18575
Expiration Date	Please enquire.