

Product Datasheet

Guinea pig Matrix Metalloproteinase 8 (MMP8) ELISA Kit (orb1204157)

Description	Guinea pig Matrix Metalloproteinase 8 (MMP8) ELISA Kit
Reactivity	Guinea pig
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	HOVPL5
Sensitivity	0.072ng/mL
Entrez	100724670
Expiration Date	Please enquire.