

Product Datasheet

Mouse Lysophosphatidylcholine Acyltransferase 4 (LPCAT4) ELISA Kit (orb1204180)

Description	Mouse Lysophosphatidylcholine Acyltransferase 4 (LPCAT4) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q6NVG1
Sensitivity	0.212ng/mL
Entrez	99010
Expiration Date	Please enquire.