

Product Datasheet

Mouse Lysyl Oxidase Like Protein 2 (LOXL2) ELISA Kit (orb1204186)

Description	Mouse Lysyl Oxidase Like Protein 2 (LOXL2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P58022
Sensitivity	0.061ng/mL
Entrez	94352
Expiration Date	Please enquire.