

Product Datasheet

Mouse IQ Motif Containing GTPase Activating Protein 2 (IQGAP2) ELISA Kit (orb1204217)

Description	Mouse IQ Motif Containing GTPase Activating Protein 2 (IQGAP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q3UQ44
Sensitivity	0.062ng/mL
Entrez	544963
Expiration Date	Please enquire.