

Product Datasheet

Guinea pig Insulin (INS) ELISA Kit (orb1204223)

Description	Guinea pig Insulin (INS) ELISA Kit
Reactivity	Guinea pig
Range	140-9000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01329
Sensitivity	64.3pg/mL
Entrez	100379579
Expiration Date	Please enquire.