

Product Datasheet

Guinea pig Interleukin 12B (IL12B) ELISA Kit (orb1204226)

Description	Guinea pig Interleukin 12B (IL12B) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q924V5
Sensitivity	6.3pg/mL
Entrez	100135554
Expiration Date	Please enquire.