

Product Datasheet

Canine Intercellular Adhesion Molecule 1 (ICAM1) ELISA Kit (orb1204231)

Description	Canine Intercellular Adhesion Molecule 1 (ICAM1) ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P33729
Sensitivity	0.064ng/mL
Entrez	403975
Expiration Date	Please enquire.