

Product Datasheet

Mouse Histone Deacetylase 2 (HDAC2) ELISA Kit (orb1204251)

Description	Mouse Histone Deacetylase 2 (HDAC2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P70288
Sensitivity	0.073ng/mL
Entrez	15182
Expiration Date	Please enquire.