

Product Datasheet

Mouse G Protein Coupled Estrogen Receptor 1 (GPER) ELISA Kit (orb1204267)

Description	Mouse G Protein Coupled Estrogen Receptor 1 (GPER) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8BMP4
Sensitivity	0.054ng/mL
Entrez	76854
Expiration Date	Please enquire.