

Product Datasheet

Mouse Glucagon Like Peptide 1 Receptor (GLP1R) ELISA Kit (orb1204274)

Description	Mouse Glucagon Like Peptide 1 Receptor (GLP1R) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O35659
Sensitivity	0.046ng/mL
Entrez	14652
Expiration Date	Please enquire.