

Product Datasheet

Mouse Growth Differentiation Factor 3 (GDF3) ELISA Kit (orb1204280)

Description	Mouse Growth Differentiation Factor 3 (GDF3) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q07104
Sensitivity	6.4pg/mL
Entrez	14562
Expiration Date	Please enquire.