

Product Datasheet

Mouse Folate Receptor 1, Adult (FOLR1) ELISA Kit (orb1204299)

Description	Mouse Folate Receptor 1, Adult (FOLR1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P35846
Sensitivity	31pg/mL
Entrez	14275
Expiration Date	Please enquire.