

Product Datasheet

Mouse Extracellular Matrix Protein 1 (ECM1) ELISA Kit (orb1204329)

Description	Mouse Extracellular Matrix Protein 1 (ECM1) ELISA Kit
Reactivity	Mouse
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q61508
Sensitivity	14.45pg/mL
Entrez	13601
Expiration Date	Please enquire.