

Product Datasheet

Mouse Dipeptidyl Peptidase 9 (DPP9) ELISA Kit (orb1204340)

Description	Mouse Dipeptidyl Peptidase 9 (DPP9) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q8BVG4
Sensitivity	25pg/mL
Entrez	224897
Expiration Date	Please enquire.