

Product Datasheet

Mouse Dickkopf Related Protein 3 (DKK3) ELISA Kit (orb1204346)

Description	Mouse Dickkopf Related Protein 3 (DKK3) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9QUN9
Sensitivity	14.3pg/mL
Entrez	50781
Expiration Date	Please enquire.