

Product Datasheet

Mouse Death Associated Protein Kinase 3 (DAPK3) ELISA Kit (orb1204358)

Description	Mouse Death Associated Protein Kinase 3 (DAPK3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	O54784
Sensitivity	0.058ng/mL
Entrez	13144
Expiration Date	Please enquire.