

Product Datasheet

Mouse Death Associated Protein Kinase 1 (DAPK1) ELISA Kit (orb1204359)

Description	Mouse Death Associated Protein Kinase 1 (DAPK1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q80YE7
Sensitivity	0.051ng/mL
Entrez	69635
Expiration Date	Please enquire.