

Product Datasheet

Mouse Cysteine Rich Protein, Angiogenic Inducer 61 (CYR61) ELISA Kit (orb1204363)

Description	Mouse Cysteine Rich Protein, Angiogenic Inducer 61 (CYR61) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P18406
Sensitivity	5.4pg/mL
Entrez	16007
Expiration Date	Please enquire.