

Product Datasheet

Mouse Cytokine Receptor Like Factor 1 (CRLF1) ELISA Kit (orb1204376)

Description	Mouse Cytokine Receptor Like Factor 1 (CRLF1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9JM58
Sensitivity	5.3pg/mL
Entrez	12931
Expiration Date	Please enquire.